

RESEARCH ARTICLE

Open Access



Using simultaneous amplification and testing method for evaluating the treatment outcome of pulmonary tuberculosis

Liping Yan[†], Heping Xiao[†] and Qing Zhang^{*}

Abstract

Background: To evaluate the utility of Simultaneous Amplification and Testing (SAT-TB) Method for monitoring anti-TB treatment response.

Methods: Serial morning sputum specimens were obtained from 377 active pulmonary tuberculosis (PTB) cases at baseline, weeks 2, months 2, 5 and 6 (newly diagnosed patients) or 8 (previously treated patients) for AmpSure assay, smear fluorescence microscopy (FM) and BACTEC MGIT 960 culture assay.

Results: After treatment of 2 weeks, sputum culture was positive in 280 patients (74.27%). Among whom, 219 patients tested positive for SAT-TB assay and 143 patients smear FM positive. The detection rate of SAT-TB (78.21%) was significantly higher than sputum FM (51.07%, $\chi^2 = 45.128$, $P < 0.001$). At the end of the second month of treatment, 157 patients (41.64%) were still culture-positive, 115 patients of them SAT-TB positive and 79 smear FM positive. The difference of detection rate between SAT-TB (73.25%) and sputum FM (50.32%) was significant ($\chi^2 = 17.480$, $P < 0.001$). When patients underwent five months of treatment, 65 patients (17.24%) with sputum culture positive was defined as treatment failure. Among whom, 60 patients (92.31%) were SAT-TB positive and 38 patients (58.46%) were smear FM positive. The detection rate of SAT-TB assay was significantly higher than sputum FM ($\chi^2 = 17.333$, $P < 0.001$).

Conclusion: Results of AmpSure assays for monitoring treatment responses can be obtained without waiting for the results of BACTEC MGIT 960 assays and most patients with treatment failures could be detected after 5 months.

Keywords: SAT-TB assay, Pulmonary tuberculosis, Culture test, Smear test, Tuberculosis treatment

Background

Tuberculosis (TB) is the leading cause of death from a single infectious disease [1]. Drug resistance further aggravates the problem. Effective drug treatment of active pulmonary tuberculosis (PTB) patients is critical for TB control. New TB patients are treated with 2HRZE/4HR regimen according to recommendation of World Health Organization (WHO) and 2HRZES/

6HRE regimen are used in previously treated patients with drug-susceptible pulmonary TB [2]. All patients in the treatment of TB should be monitored to evaluate their response to treatment by sputum examination. Monitor the patient by chest radiography is not reliable and wasteful of resources. Sputum culture conversion (SCC) is commonly used as a microbiological endpoint in the treatment of tuberculosis [3]. 2-month, 5-month and 6-month or 8-month SCC status is proxy marker for treatment outcome. However, Sputum culture is not convenient in routine clinical practice [4]. It will take weeks to obtain the results. Smear microscopy is widely used in most undeveloped

* Correspondence: zhangqingbm123@163.com

[†]Liping Yan and Heping Xiao contributed equally to this work.
Department of Tuberculosis, Shanghai Pulmonary Hospital, Tongji University School of Medicine, 507 Zhengmin Road, Shanghai 200433, People's Republic of China



countries. Sputum smear is performed at the end of the 2nd month (the intensive phase of treatment), the 5th month and treatment completion [5]. If smear is positive at month 2, sputum smear is repeated at month 3. If smear is positive at month 2 or at month 5, sputum culture and drug susceptibility testing (DST) should be obtained [5]. Smear positive at the end of the 2nd month indicates the patient may have drug-resistant *M. tuberculosis*. Smear-positivity at the end of the fifth month is defined as treatment failure and it is necessary to change the present treatment regimen. However, its sensitivity being unacceptably flawed [6] and non-viable bacteria remain visible by microscopy. In addition, microscopy could not distinguish tuberculosis mycobacteria (NTM) from *Mycobacterium tuberculosis* complex (MTBC) [5].

An increasing number of molecular methods, that is more rapid and sensitive than existing conventional tests, have become available [7–9]. Simultaneous amplification and testing methods for detection of MTBC (SAT-TB assay), that is based on nucleic acid isolation, real-time fluorescence simultaneous isothermal RNA amplification and fluorescence-labeled hybridization probes testing has the advantage of rapid, simple and good reproducibility. The results could be obtained within 120 min.

Moreover, the SAT-TB assay is intended for point-of-care testing for TB [10]. Recent researches showed that sensitivity of the SAT-TB assay for the diagnosis of PTB with sputum samples was from 39.2 to 93% [11, 12]. In patients with sputum-scarce, the sensitivity of the SAT-TB test using bronchial alveolar lavage fluid (BALF) specimens was 50.75% [13]. However, studies using the SAT-TB assay to evaluate anti-TB treatment outcome have not been reported.

For the above reasons, we designed a prospective study to investigate the efficacy of SAT-TB assays for monitoring response to anti-TB therapy with sputum specimens from new PTB patients and previously treated patients in China.

Methods

Patients

This prospective study was approved by The Ethics Committee of the Shanghai Pulmonary Hospital. Each participant gave written informed consent before enrollment.

We prospectively screened all confirmed active PTB patients based on culture positive in Shanghai Pulmonary Hospital from January 2016 to January 2017. Informations about sex; age; TB contacts; symptoms of TB; history of anti-TB treatment; comorbidities and concurrent therapies were routinely collected from each patient by attending physicians using a questionnaire. The standard 2HRZE/4HR regimen are used in newly PTB

patients according to recommendation of WHO and 2HRZES/6HRE regimen in previously treated patients under direct observed treatment (DOT) strategy [2]. The exclusion criteria included: (1) HIV test positive; (2) aged more than 17 years; (3) inability to provide sputum for examinations; (4) with concomitant extrapulmonary TB; (5) resistant to any drug of the two regimens; (6) participate in other clinical researches at the same time; (7) had a history of poor adherence in other clinical trials before.

Examinations

Morning sputum specimens were obtained from enrolled patients at baseline (at the time of diagnosis before initiation of the treatment) (T0) and at the end of 2-week (T1), 2-month (T2), 5-month (T3) and 6-month (initially diagnosed patients) or 8-month (previously treated patients) (T4) after treatment initiation. *M. tuberculosis* detection was performed using smear fluorescence microscopy (FM), bacteriological analysis and the AmpSure assay. Smear positive grade 1 was the detection threshold of a positive sputum smear. Bacteriological analysis was done using the BACTEC MGIT 960 System (Becton Dickinson Diagnostic Systems, Sparks, MD, USA), according to WHO guidelines [14]. The SAT-TB assay was performed by means of the AmpSure assay (Shanghai Rendu Biotechnology, Shanghai China) according to the protocol of manufacturer, which was described previously [15]. Initial MTB isolates obtained from patients at the time of 2-month were subjected to DST for ofloxacin (Ofx),

Isoniazid (H), rifampicin (Rfp), ethambutol (E), streptomycin (Sm), amikacin (Am), and capreomycin (Cm) using BACTEC MGIT 960 System according to WHO guidelines [14]. All of the above tests were performed in a tuberculosis reference laboratory of Shanghai Pulmonary Hospital. Quality control was performed routinely following the Manual of procedures of the National Tuberculosis Control program. Sputum specimens were given identification laboratory numbers, and technicians who carried out bacteriological investigations were blinded to the patients.

Computed tomography (CT) scans was taken on all study participants at baseline, months 2 and the end of the course or anytime if necessary. Two radiologists and a physician independently evaluated all images. A standard format was followed to make sure the radiological readings were optimized for interpretation.

Definitions

End-of-treatment outcomes were assigned according to WHO definitions [3], namely, cure: a negative sputum smear or culture result in the last month of the therapy and on previous occasion at least once; treatment

completed: completed the therapy but without a sputum smear- or culture- result in the last month of the therapy and on previous occasion at least once; failure: have a positive sputum smear or culture result at the 5th month or later or a multidrug-resistant (MDR) strain was found at anytime during the treatment; died: died due to any cause during the treatment; default: treatment was suspended for ≥ 2 consecutive months; transfer out: transferred to other registration department. Cure or treatment completion was defined as treatment success [3].

MDR is defined as resistance to at least isoniazid and rifampin. Extensively drug-resistant tuberculosis (XDR-TB), by definition, is a form of MDR-TB plus resistance to at least one drug in both of the two most significant groups of drugs in an MDR-TB schedule: fluoroquinolones (FQs) and second-line injectable agents (amikacin, capreomycin or kanamycin). Pre-XDR-TB is defined as MDR-TB plus resistance to either a FQs or a second-line injectable agent.

Statistical analyses

SPSS 18.0 software for Windows (Version 18.0, SPSS Inc., Chicago) was used to analyze the data. Numerical variables were shown as the mean \pm standard deviation. Fisher exact test or Pearson chi-squared analysis was used to analyze categorical variables. $P \leq 0.05$ was set as the level of statistical significance.

Results

We enrolled 400 suspected active PTB patients, prospectively. Twenty-three of them were excluded from our analysis because of NTM. The remaining 377 confirmed patients including 200 initially diagnosed TB cases and 177 previously treated cases with a mean age of 46.4 ± 17.5 years were eligible for this analysis and were followed during the standard anti-TB treatment (Fig. 1). These patients were predominantly males (69.5%) with a mean body mass index (BMI) of 19.0 kg/m^2 . The participants had typical CT findings that suggested active TB, including cavitory lesions, tree-in-bud

appearance and lung nodules. There are 1820 sputum specimens collected for examinations. The demographic and clinical characteristics of 377 patients were summarized in Table 1.

After treatment of 2 weeks, sputum culture was positive in 280 patients (74.27%). Among whom, 219 patients tested positive for SAT-TB assay and 143 patients smear FM positive. The detection rate of SAT-TB (78.21%) was significantly higher than sputum FM (51.07%, $X^2 = 45.128$, $P = 0.000$). At the end of the second month of treatment, 157 patients (41.64%) were still culture-positive, 115 patients of them SAT-TB positive and 79 smear FM positive. The difference of detection rate between SAT-TB (73.25%) and sputum FM (50.32%) was significant ($X^2 = 17.480$, $P = 0.000$). When patients underwent five months of treatment, 65 patients (17.24%) with sputum culture positive was defined as treatment failure. Among whom, 60 patients (92.31%) were SAT-TB positive and 40 patients (61.54%) were smear FM positive. The detection rate of SAT-TB assay was significantly higher than sputum FM ($X^2 = 17.333$, $P < 0.001$). Ultimately, 312 patients (82.76%) cured with culture clearing and imaging improvement. The results of SAT-TB assay, smear FM, and BACTEC MGIT 960 culture at different time points were summarized in Table 2.

The overall sensitivity of SAT-TB assay was 86.92%, which was significantly higher than that of smear FM (56.43%) ($P < 0.001$). The overall specificity of SAT-TB assay was 100%, which was also higher than that of smear FM (99.69%) ($P = 0.156$). The PPV and NPV of SAT-TB assay was 100% and 84.55%, respectively. The PPV and NPV of smear FM was 99.6% and 62.28%, respectively. The difference between PPV had no statistical significance ($P = 0.080$), but the difference between NPV had statistical significance ($P < 0.001$) (Fig. 2).

Initial MTB isolates obtained from positive culture of 157 patients at the end of the second month of treatment were subjected to DST. The DST results were summarized in Table 3. There were 6 XDR-TB, 9 pre-XDR-TB, 32 MDR-TB and 6 mono-drug resistant (resistant to




Table 1 Patient characteristics

| Variable | New PTB (n = 200) | Previously treated (n = 177) | Total (n = 377) |
|----------------------------------|----------------------|---------------------------------|--------------------|
| Men (n [%]) | 135 (67.5) | 127 (71.75) | 262 |
| Women (n [%]) | 65 (32.5) | 50 (28.25) | 115 |
| Age (y) (median) (range) | 44.6 ± 17.7 | 47.4 ± 19.7 | 46.4 ± 17.5 |
| Body mass index (median) (range) | 19.1 (13–27) | 18.9 (13–27) | |
| Cured cases | 173 | 139 | 312 |
| Failed cases | 27 | 38 | 65 |

streptomycin) TB patients. All of the 47 MDR/XDR-TB and 3 mono-drug resistant as well as 15 susceptible TB patients failed.

Discussion

In the current study, we used SAT-TB assay to evaluate treatment outcome in Chinese PTB patients compared with smear FM and BACTEC MGIT 960 at different time points, and observed that the SAT-TB assay is preferred. All the patients that was in the treatment of TB should be monitored the response to treatment. If sputum smear is positive at month 2, sputum smear is repeated at month 3. If smear is positive at month 3, sputum culture and DST should be obtained. Smear-positivity or culture-positivity at the end of the fifth month is defined as treatment failure [3]. Unfortunately, sputum smear and culture have significant limitations. The technological development of nucleic acid amplification has led to advances in the early detection of tubercle bacilli. The GeneXpert assay (Cepheid, Sunnyvale, CA) was introduced currently in China. The high price of the GeneXpert assay and the requirement of costly especial instrument for the MTD test (Gen-Probe; San Diego, CA, USA) prevent these assays from being a routine test in resource-limited countries [16,

17]. The SAT-TB assay has another advantage over PCR, as the detection objective is rRNA. RNA is much more unstable than DNA, so a positive result indicates the presence of viable MTBC. Therefore the false positive rate decreased. In addition, the SAT-TB assay can be performed on real-time PCR instruments, which can be found in most clinical laboratories.

The SAT-TB assay was demonstrated to have better sensitivity and high specificity for the detection of tubercle bacilli in our study. Over the intensive phase of therapy, 73.25% of culture-positive patients could be detected with the SAT-TB assay within 2 h. Month-2 culture conversion is generally used as substitute index of anti-TB treatment efficiency [18]. Detection of tubercle bacilli at the end of the intensive phase remains important as it may indicate that the patient may have multidrug-resistant *M. tuberculosis* and DST should be performed. If we use sputum smear for monitoring the treatment response, only about half of culture-positive patients (50.32%) could be detected. If we use the SAT-TB assay, DST could be performed in three-quarters of these patients. The number of patients in need of DST was increased by about a quarter. Drug-resistant tuberculosis (DR-TB) should be identified as early as possible not only for achieving the best treatment outcome but also for controlling the spread of DR-TB. Smear-positivity at the end of the fifth month is defined as treatment failure and it is necessary to change the present treatment regimen. In this study, 65 patients with sputum culture positive were defined as treatment failure. Among whom, 60 patients were SAT-TB assay positive. This means that if we use the SAT-TB assay for monitoring the treatment response, 92.31% of patients with treatment failure can be found timely and changed to appropriate therapy. However, only 58.46% of these culture-positive patients were sputum smear positive. After a month or three, physicians may mistakenly declare the treatment successful for 41.54% of these failure patients. It's worth mentioning that two sputum smear

Table 2 Summary of SAT-TB, smear FM, and culture results

| | | SAT-TB (+) | | SAT-TB (-) | |
|------------------|----------------|--------------|--------------|------------|--------------|
| | | Smear (+) | Smear (-) | Smear (+) | Smear (-) |
| before treatment | Mt-culture (+) | 236 (62.6%) | 134 (35.55%) | 0 | 7 (1.86%) |
| | Mt-culture (-) | 0 | 0 | 0 | 0 |
| 2nd week | Mt-culture (+) | 143 (51.07%) | 76 (27.15%) | 0 | 61 (21.78%) |
| | Mt-culture (-) | 0 | 0 | 0 | 97 (25.73%) |
| 2nd month | Mt-culture (+) | 79 (50.32%) | 36 (22.93%) | 0 | 42 (26.75%) |
| | Mt-culture (-) | 0 | 0 | 0 | 220 (59.46%) |
| 5th month | Mt-culture (+) | 38 (58.46%) | 22 (33.84%) | 0 | 5 (7.7%) |
| | Mt-culture (-) | 0 | 0 | 2 (0.65%) | 310 (99.35%) |
| 6th/8th month | Mt-culture (+) | 0 | 0 | 0 | 0 |
| | Mt-culture (-) | 0 | 0 | 0 | 312 (82.76%) |




Fig. 2 ROC curve for SAT-TB and smear FM

positive were identified as false-positives or laboratory contamination according to culture and SAT-TB assay negative and imaging improvement at the end of the fifth month. Ultimately, the treatment success rate in this study was 82.76%, which was comparable to the average success rate in the world. Globally, the treatment success rate for the 5.9 million new and relapse cases was 83% in 2015 [1]. False-negative findings, however, are inevitable with the current technology [19–21]. In addition, the major barrier in applying molecular detection method in evaluating treatment response is the problem of false positive rather than false negative. In our research, there was no false positive (sputum culture-negative but SAT-TB assay positive) for SAT-TB assay.

Treatment outcome of anti-TB treatment in bacteriologically confirmed PTB patients is evaluated principally by serial bacteriologic examinations [2], whereas Treatment outcome in clinically diagnosed PTB patients are usually monitored clinically or radiographically. To overcome these

shortcomings, numerous biomarkers have been evaluated [4, 22], although most of the markers do not have high levels of validity. More and more studies confirmed that the possible function of plasma-based IFN- γ -release assays (IGRAs) for evaluation of reaction to anti-TB therapy is objectionable [23]. The IGRA is less useful in terms of treatment monitoring, given that the IGRA is based on the TB antigen-stimulated immunologic reaction, which tends to persist and shows blunted fluctuations during the longitudinal monitoring of TB [24]. Lago et al. observed a consistent decrease in IL-10 levels in active PTB patients at all times of therapy, suggesting that this cytokine could be used as a helpful biomarker for evaluating disease progression [25]. Therefore, novel accessing means of anti-TB therapy outcome are required to make anti-TB therapy better in some clinical circumstances.

In this study, 65% of patients were male. Men seem to be more affected with TB than women with a male/female ratio of 1.86 for the worldwide case notification rate in 2017 [26]. This excess of male pulmonary TB cases is seen in almost all countries of the world, and this gender effect is thought to be related to many factors, such as access to healthcare, nutritional status, socioeconomic and cultural factors. TB disease is not gender specific. Male and female are equally susceptible.

There are several limitations in this study. First, the number of patients with microbiological diagnosed active TB is relatively small. Second, this study was limited to only HIV-negative PTB patient group, which limits its generalization to the HIV-infected group.

Table 3 Summary of drug susceptibility testing results of 157 patients

| | | XDR | Pre-XDR | MDR | Sm-resistant | Susceptible | Total |
|--------------------|---------|-----|---------|-----|--------------|-------------|-------|
| New PTB | cure | 0 | 0 | 0 | 1 | 34 | 35 |
| | failure | 0 | 2 | 9 | 1 | 15 | 27 |
| Previously treated | cure | 0 | 0 | 0 | 2 | 55 | 57 |
| | failure | 6 | 7 | 23 | 2 | 0 | 38 |
| Total | | 6 | 9 | 32 | 6 | 104 | 157 |

HIV-positive cases, a clinically significant population, were removed from this study as they must be shifted to special hospitals in China. Third, a few variables such as chest imaging manifestations were not explored in this analysis. Finally, culture after treatment was not followed up in this study. Despite these shortcomings, as far as we know, this is likely to be the first prospective study to evaluate the utility of SAT-TB assays to access TB therapy outcome.

Conclusions

In summary, our research suggests that the SAT-TB assay is preferred to evaluate treatment outcome in Chinese PTB patients compared with smear FM and BACTEC MGIT 960 at different time points. Results of AmpSure assays for monitoring treatment responses can be obtained without waiting for the results of BACTEC MGIT 960 assays and most patients with treatment failures could be detected after 5 months. Whether SAT-TB assays can be used to access TB therapy outcome needs to be explored further.

Abbreviations

Am: Amikacin; BALF: Bronchial alveolar lavage fluid; Cm: Capreomycin; CT: Computed tomography; DOT: Direct observed treatment; DR-TB: Drug-resistant tuberculosis; DST: Drug susceptibility testing; E: Ethambutol; FM: Fluorescence microscopy; FQs: Fluoroquinolones; H: Isoniazid; IGRAs: IFN- γ -release assays; MDR: Multidrug-resistant; MTBC: Mycobacterium tuberculosis complex; NTM: Non-tuberculous mycobacteria; Ofx: Ofloxacin; PTB: Pulmonary tuberculosis; Rfp: Rifampicin; SAT-TB assay: Simultaneous amplification and testing methods for detection of MTBC; SCC: Sputum culture conversion; Sm: Streptomycin; TB: Tuberculosis; WHO: World Health Organization; XDR-TB: Extensively drug-resistant tuberculosis

Funding

This work was supported by a grant from the National Key Scientific and Technical Program (Grant No. 2013ZX10003009).

Availability of data and materials

All data generated or analyzed during this study are included in this published article. The datasets used and/or analysed during the current study are available from the corresponding author on reasonable request.

Authors' contributions

QZ and HX were responsible for the conception and design of the study. All author were responsible for acquisition and analysis of data; furthermore, LY and QZ were in charge of statistical analysis. All authors drafted the manuscript; QZ revised and approved the final version of the manuscript. All authors read and approved the final manuscript.

Ethics approval and consent to participate

This prospective study was approved by The Ethics Committee of the Shanghai Pulmonary Hospital. Each participant gave written informed consent before enrollment.

Consent for publication

Not applicable.

Competing interests

The authors declare that they have no competing interests.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Received: 7 May 2018 Accepted: 27 September 2018

Published online: 11 October 2018

References

1. WHO: Global Tuberculosis Report 2017. In Geneva: World Health Organization; 2017.
2. WHO. Guidelines for treatment of drug-susceptible tuberculosis and patient care, 2017 update. Geneva: World Health Organization; 2017.
3. Wallis RS, Kim P, Cole S, Hanna D, Andrade BB, Maeurer M, et al. Tuberculosis biomarkers discovery: developments, needs, and challenges. *Lancet Infect Dis*. 2013;13(4):362–72.
4. Ritchie SR, Harrison AC, Vaughan RH, Calder L, Morris AJ. New recommendations for duration of respiratory isolation based on time to detect *Mycobacterium tuberculosis* in liquid culture. *Eur Respir J*. 2007;30(3):501–7.
5. Eurosurveillance editorial team. WHO revised definitions and reporting framework for tuberculosis. *Euro Surveill*. 2013;18(16):20455.
6. Desikan P. Sputum smear microscopy in tuberculosis: is it still relevant? *Indian J Med Res*. 2013;137(3):442–4.
7. Dinnes J, Deeks J, Kunst H, Gibson A, Cummins E, Waugh N, et al. A systematic review of rapid diagnostic tests for the detection of tuberculosis infection. *Health Technol Assess*. 2007;11(3):1–196.
8. Piersimoni C, Scarparo C, Piccoli P, Rigon A, Ruggiero G, Nista D, et al. Performance assessment of two commercial amplification assays for direct detection of *Mycobacterium tuberculosis* complex from respiratory and extrapulmonary specimens. *J Clin Microbiol*. 2002;40(11):4138–42.
9. Yan L, Xiao H, Zhang Q. Systematic review: comparison of Xpert MTB/RIF, LAMP and SAT methods for the diagnosis of pulmonary tuberculosis. *Tuberculosis*. 2016;96:75–86.
10. Mdivani N, Li H, Akhalaia M, Gegia M, Goginashvili L, Kernodle DS, et al. Monitoring therapeutic efficacy by real-time detection of *Mycobacterium tuberculosis* mRNA in sputum. *Clin Chem*. 2009;55(9):1694–700.
11. Cui Z, Wang Y, Fang L, Zheng R, Huang X, Liu X, et al. Novel real-time simultaneous amplification and testing method to accurately and rapidly detect *Mycobacterium tuberculosis* complex. *J Clin Microbiol*. 2012;50(3):646–50.
12. Fan L, Zhang Q, Cheng L, Liu Z, Ji X, Cui Z, et al. Clinical diagnostic performance of the simultaneous amplification and testing methods for detection of the *Mycobacterium tuberculosis* complex for smear-negative or sputum-scarce pulmonary tuberculosis in China. *Chin Med J*. 2014;127(10):1863–7.
13. Yan L, Zhang Q, Xiao H. Clinical diagnostic value of simultaneous amplification and testing for the diagnosis of sputum-scarce pulmonary tuberculosis. *BMC Infect Dis*. 2017;17(1):545.
14. WHO. Guidelines for surveillance of drug resistance in tuberculosis. 4th ed. Geneva: World Health Organization; 2009.
15. Yan L, Tang S, Yang Y, Shi X, Ge Y, Sun W, et al. A large cohort study on the clinical value of simultaneous amplification and testing for the diagnosis of pulmonary tuberculosis. *Medicine*. 2016;95(4):e2597.
16. Du J, Huang Z, Luo Q, Xiong G, Xu X, Li W, et al. Rapid diagnosis of pleural tuberculosis by Xpert MTB/RIF assay using pleural biopsy and pleural fluid specimens. *J Res Med Sci*. 2015;20(1):26–31.
17. Jonas V, Alden MJ, Curry JI, Kamisango K, Knott CA, Lankford R, et al. Detection and identification of *Mycobacterium tuberculosis* directly from sputum sediments by amplification of rRNA. *J Clin Microbiol*. 1993;31(9):2410–6.
18. Mitchison D, Davies G. The chemotherapy of tuberculosis: past, present and future. *Int J Tuberc Lung Dis*. 2012;16(6):724–32.
19. Boehme CC, Nabeta P, Hillemann D, Nicol MP, Shenai S, Krapp F, et al. Rapid molecular detection of tuberculosis and rifampin resistance. *N Engl J Med*. 2010;363(11):1005–15.
20. Boehme CC, Nicol MP, Nabeta P, Michael JS, Gotuzzo E, Tahirli R, et al. Feasibility, diagnostic accuracy, and effectiveness of decentralised use of the Xpert MTB/RIF test for diagnosis of tuberculosis and multidrug resistance: a multicentre implementation study. *Lancet*. 2011;377(9776):1495–505.
21. Sethi S, Singh S, Dhatwalia SK, Yadav R, Mewara A, Singh M, et al. Evaluation of in-house loop-mediated isothermal amplification (LAMP) assay for rapid diagnosis of *M. tuberculosis* in pulmonary specimens. *J Clin Lab Anal*. 2013;27(4):272–6.
22. Mihret A, Bekele Y, Bobosha K, Kidd M, Aseffa A, Howe R, et al. Plasma cytokines and chemokines differentiate between active disease and non-active tuberculosis infection. *J Infect*. 2013;66(4):357–65.

23. Chiappini E, Fossi F, Bonsignori F, Sollai S, Galli L, de Martino M. Utility of interferongamma release assay results to monitor anti-tubercular treatment in adults and children. *Clin Ther.* 2012;34(5):1041–8.
24. Denkinger CM, Pai M, Patel M, Menzies D. Gamma interferon release assay for monitoring of treatment response for active tuberculosis: an explosion in the spaghetti factory. *J Clin Microbiol.* 2013;51(2):607–10.
25. Lago PM, Boéchat N, Migueis DP, Almeida AS, Lazzarini LC, Saldanha MM, et al. Interleukin-10 and interferon-gamma patterns during tuberculosis treatment: possible association with recurrence. *Int J Tuberc Lung Dis.* 2012; 16(5):656–9.
26. Khan KB. Understanding the gender aspects of tuberculosis: a narrative analysis of the lived experiences of women with TB in slums of Delhi, India. *Health Care Women Int.* 2012;33(1):3–18.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

